

PRODUCT-SPECIFIC PROTOCOLS

WESTERN BLOT (CL488-15105)

Sample type	Amount of protein loaded	Membrane	Transfer type	Blocking buffer	Primary antibody dilution	Incubation time	Detection method
HeLa cells	50ug	PVDF	Wet	5% milk in TBST	1:500	37°C 1 hour	Image System

PROTOCOL

1. Prepare sample lysate, heat lysate in sample buffer at 100°C for 5min, and resolve proteins via SDS-PAGE. (Note: Prepare sample lysate without heating or heating at 37°C according to picture recommended for some membrane proteins.)
2. Transfer proteins from the gel onto the membrane.
3. Incubate membrane with Blocking buffer on a rocking platform.
4. Prepare the primary antibody in Blocking buffer.
5. Incubate membrane with primary antibody on a rocking platform.
6. Wash the membrane 3 times for 10 minutes each in 1X TBST
7. Expose the membrane in fluorescence image system, such as iBright. The exposure time can be optimized.

Visit ptglab.com for full product range