

Influenza A H3N2 (A/Brisbane/10/2007) Hemagglutinin / HA Antibody, Rabbit MAb

Catalog Number: 11056-R006



Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Influenza A H3N2 (A/Brisbane/10/2007) Hemagglutinin / HA Protein (Catalog#11056-V08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Influenza A H3N2 (A/Brisbane/10/2007) Hemagglutinin / HA (Catalog#11056-V08H; ABW23353.1; Met1-Ile531).
Ig Type:	Rabbit IgG
Clone ID:	006
Specificity:	Influenza A H3N2 (A/Brisbane/10/2007) Hemagglutinin / HA
	Has cross-reactivity in ELISA with H3N2 (A/Darwin/9/2021) HA Protein (Cat# 40859-V08B) H3N2 (A/Darwin/9/2021) HA Protein (Cat# 40859-V08H) H3N2 (A/Darwin/1415/2025) HA trimer Protein (Cat# 41128-V08H) H3N2 (A/Massachusetts/18/2022; A/Thailand/8/2022) HA Protein (Cat# 40992-V08B) H3N2 (A/Massachusetts/18/2022; A/Thailand/8/2022) HA Protein (Cat# 40992-V08H) No cross-reactivity in ELISA with H1N1 (A/Wisconsin/67/2022) HA Trimer Protein (Cat# 40940-V08H2) H1N1 (A/Wisconsin/67/2022) HA Protein (Cat# 40940-V08B) H1N1 (A/Wisconsin/67/2022) HA Protein (Cat# 40940-V08H) Influenza B (B/Austria/1359417/2021) HA Protein (Cat# 40862-V08H) Influenza B (B/Austria/1359417/2021) HA Protein (Cat# 40862-V08B)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2 °C -8 °C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20 °C to -80 °C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	Hemagglutinin,HA
APPLICATIONS	
Applications:	ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Influenza A H3N2 (A/Brisbane/10/2007) Hemagglutinin / HA.

Please Note: Optimal concentrations/dilutions should be determined by the end user.