

APC-Labeled Monoclonal Anti-Human CD2 Antibody, Mouse IgG1 (RPA-2.10) (0.03% Proclin)

Catalog # CD2-AHFM8



Source

Anti-Human CD2 Antibody, Mouse IgG1 (RPA-2.10) is a recombinant mouse monoclonal antibody expressed in HEK293 cells, engineered with a humanized segment in its Fab heavy chain, which provides higher batch consistency and long term security of supply.

Clone

RPA-2.10

Isotype

Mouse IgG1, Kappa

Host Species

Mouse

Reactivity

Human

Immunogen

Purified Human CD2 Protein

Conjugate

APC

Excitation Wavelength: 640 nm

Emission Wavelength: 661 nm

Specificity

This product is a specific antibody specifically reacts with CD19 protein.

Application

Flow Cytometry (Evaluation of the expression of CD2 on Human cells).

Recommended Dilution

1:50

Form

Lyophilized

Formulation

Lyophilized from a 0.22 µm-filtered solution in PBS (pH 7.4) containing 0.03% ProClin 300 and 0.2% BSA, with trehalose as protectant.

Please contact us for customized product forms or formulations.

Reconstitution

Please refer to the Certificate of Analysis (CoA) for specific instructions.

Shipping

Lyophilized product is shipped at ambient temperature.

Storage

For long term storage, the product should be stored in a lyophilized state at -20°C or lower.

Please protect from light and avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 24 months in lyophilized state;
- -70°C for 12 months after reconstitution.
- 2-8 °C for 12 months after reconstitution.


Notices

Product Specific Notices: For research use only.

ACRO Quality Management System

- [QMS\(ISO, GMP\)](#)
- [Quality Advantages](#)
- [Quality Control Process](#)

Bioactivity-FACS



Flow cytometric analysis of Human peripheral blood lymphocytes staining with APC-Labeled Monoclonal Anti-Human CD2 Antibody, Mouse IgG1 (RPA-2.10) (Cat. No. CD2-AHFM8) at 1:50 dilution (2 μ L of the antibody stock solution corresponds to labeling of 1e6 PBMCs in a final volume of 100 μ L), compared with APC Mouse IgG1, κ Isotype Ctrl Antibody. APC signal was used to evaluate the binding activity (QC tested).

Background

T-cell surface antigen CD2 is also known as Erythrocyte receptor, LFA-2, LFA-3 receptor, Rosette receptor, T-cell surface antigen T11/Leu-5 and SRBC, is a single-pass type I membrane protein found on the surface of T cells and natural killer (NK) cells. CD2 is a member of the immunoglobulin superfamily. CD2 / SRBC contains 1 Ig-like C2-type (immunoglobulin-like) domain and 1 Ig-like V-type (immunoglobulin-like) domain. CD2 / SRBC interacts with other adhesion molecules, such as lymphocyte function-associated antigen-3 (LFA-3 / CD58) in humans, or CD48 in rodents, which are expressed on the surfaces of other cells. In addition to its adhesive properties, CD2 also acts as a co-stimulatory molecule on T and NK cells. CD2 is a specific marker for T cells and NK cells, and can therefore be used in immunohistochemistry to identify the presence of such cells in tissue sections.

Discounts, Gifts,
and more!



www.acrobiosystems.com