




Recombinant Human C-C motif chemokine 7 protein (CCL7) (Active)

Product Code	CSB-AP000881HU
Abbreviation	Recombinant Human CCL7 protein (Active)
Uniprot No.	P80098
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered 20 mM PB, pH 7.4, 150 mM NaCl
Product Type	Chemokine
Immunogen Species	Homo sapiens (Human)
Biological Activity	Fully biologically active when compared to standard. The biologically active determined by a chemotaxis bioassay using human monocytes is in a concentration range of 10-100 ng/ml.
Purity	≥97% as determined by SDS-PAGE.
Sequence	QPVGINTSTT CCYRFINKKI PKQRLESYRR TTSSHCPREA VIFKTKLDKE ICADPTQKWV QDFMKHLDKK TQTPKL
Research Area	Immunology
Source	E.coli
Target Names	CCL7
Expression Region	24-99aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	9.0 kDa
Protein Length	Full Length of Mature Protein
PubMed ID	8461011; 7916328; 8318676; 15489334; 1613466; 15340161; 8898111; 9109648;

Image





Endotoxin

Less than 1.0 EU/μg as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.