



FOS Monoclonal Antibody

Product Code	CSB-MA080264
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01100
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgG
Alias	c-Fos
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOS
Clone No.	13H4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.