



TAS2R45 Antibody

Product Code	CSB-PA008596
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59539
Immunogen	Synthesized peptide derived from the C-terminal region of Human T2R45.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TAS2R45; GPR59; Taste receptor type 2 member 45; T2R45; G-protein coupled receptor 59
Immunogen Species	Homo sapiens (Human)
Target Names	TAS2R45
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.