



GABBR2 Antibody

Product Code	CSB-PA009135ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75899
Immunogen	Recombinant Human Gamma-aminobutyric acid type B receptor subunit 2 protein (60-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-aminobutyric acid type B receptor subunit 2 (GABA-B receptor 2) (GABA-B-R2) (GABA-BR2) (GABABR2) (Gb2) (G-protein coupled receptor 51) (HG20), GABBR2, GPR51 GPRC3B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GABBR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.