



pInA Antibody

Product Code	CSB-PA302126ZA01LMS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80214
Immunogen	Recombinant Lactobacillus plantarum pInA protein
Raised In	Rabbit
Species Reactivity	Lactobacillus plantarum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Lactobacillus plantarum (strain ATCC BAA-793 / NCIMB 8826 / WCFS1)
Lead Time	Made-to-order (12-14 weeks)
Target Names	pInA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.