



2a Antibody

Product Code	CSB-PA325583ZA01BJN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18517
Immunogen	Recombinant Bovine coronavirus 2a protein
Raised In	Rabbit
Species Reactivity	Bovine coronavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Bovine coronavirus (strain Quebec) (BCoV) (BCV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	2a
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.