



Thanatin Antibody

Product Code	CSB-PA345509ZA01PNT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55788
Immunogen	Recombinant Podisus maculiventris Thanatin protein
Raised In	Rabbit
Species Reactivity	Podisus maculiventris
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Podisus maculiventris (Spined soldier bug)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Thanatin
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.