



HA Antibody

Product Code	CSB-PA356317LA01IDG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04664
Immunogen	Recombinant Influenza A virus Hemagglutinin protein (1-328AA)
Raised In	Rabbit
Species Reactivity	Influenza A virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemagglutinin [Cleaved into: Hemagglutinin HA1 chain] (Fragment), HA
Immunogen Species	Influenza A virus (strain A/England/878/1969 H3N2)
Research Area	Others
Target Names	HA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.