



PRO2 Antibody

Product Code	CSB-PA521315ZA01EUQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O24650
Immunogen	Recombinant Phleum pratense PRO2 protein
Raised In	Rabbit
Species Reactivity	Phleum pratense
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Phleum pratense (Common timothy)
Lead Time	Made-to-order (12-14 weeks)
Target Names	PRO2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.