



N Antibody

Product Code	CSB-PA700347ZA01HIY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4VID0
Immunogen	Recombinant Human coronavirus OC43 N protein
Raised In	Rabbit
Species Reactivity	Human coronavirus OC43
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Human coronavirus OC43 (HCoV-OC43)
Lead Time	Made-to-order (12-14 weeks)
Target Names	N
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.