



N Antibody

Product Code	CSB-PA889516ZA01BJI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9QAR1
Immunogen	Recombinant Bovine coronavirus N protein
Raised In	Rabbit
Species Reactivity	Bovine coronavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Bovine coronavirus (strain LSU-94LSS-051) (BCoV-LSU) (BCV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	N
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.