




Eribulin Recombinant Monoclonal Antibody


Product Code	CSB-RA943989I1m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Eribulin-BSA
Species Reactivity	Chemical
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG2a
Clonality	Monoclonal
Product Type	Recombinant Antibody
Purity	≥ 95% as determined by SEC-HPLC.
Research Area	Others
Target Names	Eribulin
Clone No.	45C2G9

Image



The Binding Activity of Farletuzumab ecteribulin (Synonyms: MORAb-202) with Anti-Eribulin antibody

Activity: Measured by its binding ability in a functional ELISA. Immobilized Farletuzumab ecteribulin (Synonyms: MORAb-202) at 2 µg/mL can bind Anti-Eribulin antibody, the EC₅₀ is 0.8548 to 1.015 ng/mL.



The purity of Eribulin was greater than 95% as determined by SEC-HPLC

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.