


Biotinylated Human LILRA3 (C-6His-Avi)(Cat.No.:CY50)

Source:	Human Cells
Description:	Biotinylated Recombinant Human Leukocyte Immunoglobulin-Like Receptor Subfamily A Member 3 is produced by our Mammalian expression system and the target gene encoding Gly24-Glu439 is expressed with a 6His, Avi tag at the C-terminus.
Accession:	AAH28208.1
Known As:	Leukocyte immunoglobulin-like receptor subfamily A member 3; CD85 antigen-like family member E; Immunoglobulin-like transcript 6; ILT-6; Leukocyte immunoglobulin-like receptor 4; LIR-4 and Monocyte inhibitory receptor HM43/HM31
Predicted Mol Mass:	47.9 KDa
Apparent Mol Mass:	70-90 KDa, reducing conditions
Endotoxin:	< 1 EU/μg as determined by LAL test.
Formulation:	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Reconstitution:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Shipping:	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Storage:	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Background:	Leukocyte immunoglobulin-like receptor subfamily A member 3 is also known as CD85 antigen-like family member E, Immunoglobulin-like transcript 6, ILT-6, Leukocyte immunoglobulin-like receptor 4, LIR-4 and Monocyte inhibitory receptor HM43/HM31. In humans, it is encoded by the LILRA3 gene. It acts as soluble receptor for class I MHC antigens. Binds both classical and non-classical HLA class I molecules but with reduced affinities compared to LILRB1 or LILRB2. It is detected in B-cells, natural killer (NK) cells, peripheral blood monocytes and lung.

Version 26.2.3

Purity-SDS-PAGE:



Greater than 95% as determined by reducing SDS-PAGE. (QC verified)