

G-CSF/CSF3 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10007-MM32

GENERAL INFORMATION

Immunogen:	Recombinant Human G-CSF/CSF3 Protein (Catalog#10007-HNCE)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human G-CSF/CSF3 (rh G-CSF/CSF3; Catalog#10007-HNCE; NP_757373.1; Ala30-Pro204). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	32
Specificity:	Human G-CSF/CSF3
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human G-CSF/CSF3.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.