

KRT14 / Cytokeratin 14 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 100129-RP02

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to the C-terminus of the human Cytokeratin 14.
Preparation	Produced in rabbits immunized with purified, a synthetic peptide corresponding to the C-terminus of the human Cytokeratin 14, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human, Mouse, Rat
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB, IHC-P
---------------	-----------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:500-1:2000
Western Blot	WB: 1:500-1:2000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

KRT14 / Cytokeratin 14 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist


Catalog Number: 100129-RP02




Immunohistochemical staining CK14 in human skin with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).




Immunohistochemical staining CK14 in human esophageal carcinoma with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).



Immunohistochemical staining CK14 in mouse skin with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).



Immunohistochemical staining CK14 in rat skin with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).



Anti-Cytokeratin14 rabbit polyclonal antibody at 1:500 dilution
Lane A: A431 Whole Cell Lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG H&L (Dylight800) at 1/10000 dilution.

Developed using the Odyssey technique.
Performed under reducing conditions.

Predicted band size: 52 kDa

Observed band size: 52 kDa (We are unsure as to the identity of these extra bands.)