



Catalog Number: 100140-R012

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to the center region of the human GFAP
Preparation	This antibody was obtained from a rabbit immunized with a synthetic peptide corresponding to the center region of the human GFAP.
Ig Type:	Rabbit IgG
Clone ID:	012
Specificity:	Human GFAP
Formulation:	0.2 µm filtered solution in PBS with 10% Trehalose, 0.02% Tween 80, pH7.0
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB
----------------------	----

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:500-1:1000
---------------------	------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.


GFAP Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 100140-R012

EliteRmab® is a registered trademark of Sino Biological Inc.



Anti-GFAP rabbit monoclonal antibody at
1:500 dilution

Lane A: Mouse brain tissue lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG H&L (Dylight800) at
1/10000 dilution.

Developed using the Odyssey technique.
Performed under reducing conditions.

Predicted band size: 50 kDa

Observed band size: 45 kDa