

ARID3B Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 100355-T32

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to the N-terminus of the human ARID3B
Preparation	Produced in rabbits immunized with A synthetic peptide corresponding to the N-terminus of the human ARID3B , and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human ARID3B
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB
---------------	----

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:500-1:1000
--------------	------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

ARID3B Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 100355-T32



Anti-ARID3B rabbit polyclonal antibody at
1:500 dilution
Lane A: 293T Whole Cell Lysate

Lysates/proteins at 30 µg per lane.
Secondary
Goat Anti-Rabbit IgG H&L (Dylight800) at
1/10000 dilution.

Developed using the Odyssey technique.
Performed under reducing conditions.

Predicted band size: 61 kDa
Observed band size: 61 kDa
(We are unsure as to the identity of these
extra bands.)