

CD27 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10039-T24

GENERAL INFORMATION	
Immunogen:	Recombinant Human CD27 protein (Catalog#10039-H31H)
Preparation	Produced in rabbits immunized with purified, recombinant Human CD27 (rh CD27; Catalog#10039-H31H; P26842; Met1-Ile192). CD27 specific IgG was purified by Human CD27 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CD27
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA, IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:1000-1:4000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human CD27.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD27 Antibody, Rabbit PAb, Antigen Affinity Purified


Catalog Number: 10039-T24



Sino Biological
Biological Solution Specialist



Immunochemical staining of human CD27 in human spleen with rabbit polyclonal antibody (1:2000, formalin-fixed paraffin embedded sections).



Immunochemical staining of human CD27 in human tonsil with rabbit polyclonal antibody (1:2000, formalin-fixed paraffin embedded sections).