

CD27/TNFRSF7 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10039-T56

GENERAL INFORMATION	
Immunogen:	Recombinant Human CD27/TNFRSF7 Protein (Catalog#10039-H31H)
Preparation	Produced in rabbits immunized with purified, recombinant Human CD27/TNFRSF7 (rh CD27/TNFRSF7; Catalog#10039-H31H; P26842; Met1-Ile192). CD27/TNFRSF7 specific IgG was purified by Human CD27/TNFRSF7 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CD27/TNFRSF7
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	WB, ELISA, IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:2000-1:10000
Western Blot	WB: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD27/TNFRSF7 Antibody, Rabbit PAb, Antigen Affinity Purified


Catalog Number: 10039-T56



Sino Biological
Biological Solution Specialist



Immunohistochemical staining of human CD27 in human spleen with rabbit polyclonal antibody at 1:5000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human CD27 in human tonsil with rabbit polyclonal antibody at 1:5000 dilution, formalin-fixed paraffin embedded sections.



Anti-CD27 rabbit polyclonal antibody at 1:1000 dilution

Lane A: Raji Whole Cell Lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 29 kDa

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>