

CD13/ANPEP Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10051-R422

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Human CD13/ANPEP Protein (Catalog#10051-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CD13/ANPEP (rh CD13/ANPEP; Catalog#10051-H08H; NP_001141.2; Lys69-Lys967,603Met/Ile).
Ig Type:	Rabbit IgG
Clone ID:	422
Specificity:	Human CD13/ANPEP
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	IHC-P
---------------	-------

RECOMMENDED CONCENTRATION


IHC-P	IHC-P: 1:100-1:500
-------	--------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.




Catalog Number: 10051-R422

EliteRmab® is a registered trademark of Sino Biological Inc.



Immunohistochemical staining of human ANPEP in human kidney with rabbit monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human ANPEP in human hepatoma with rabbit monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.