

CAV3 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 100547-T08

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to the N-terminus of the Human CAV3
Preparation	Produced in rabbits immunized with A synthetic peptide corresponding to the N-terminus of the Human CAV3, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human, Mouse, Rat
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	IHC-P
---------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:1000-1:5000
-------	----------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CAV3 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist


Catalog Number: 100547-T08




Immunohistochemical staining CAV3 in human skeletal muscle with rabbit polyclonal antibody (1:2000, formalin-fixed paraffin embedded sections).




Immunohistochemical staining CAV3 in human heart with rabbit polyclonal antibody (1:2000, formalin-fixed paraffin embedded sections)




Immunohistochemical staining CAV3 in mouse heart with rabbit polyclonal antibody at 1:2000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining CAV3 in mouse skeletal muscle with rabbit polyclonal antibody at 1:2000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining CAV3 in rat heart with rabbit polyclonal antibody at 1:2000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining CAV3 in rat skeletal muscle with rabbit polyclonal antibody at 1:2000 dilution, formalin-fixed paraffin embedded sections.