

DPP6 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 100592-T36

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to the center region of the Human DPP6
Preparation	Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human DPP6, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human DPP6
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB, IP
---------------	--------

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:500-1:2000
IP	IP:1-5uL/mg of lysate


Please Note: Optimal concentrations/dilutions should be determined by the end user.

DPP6 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 100592-T36



DPP6 was immunoprecipitated using:
Lane A: 0.5 mg Rat brain Whole Cell Lysate


4 μ L anti-DPP6 rabbit polyclonal antibody
and 60 μ g of Immunomagnetic beads
Protein A/G.

Primary antibody:
Anti-DPP6 rabbit polyclonal antibody, at
1:100 dilution

Secondary antibody:
Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000
dilution

Developed using the ECL technique.
Performed under reducing conditions.

Predicted band size: 98 kDa
Observed band size: 115 kDa



Anti-DPP6 rabbit monoclonal antibody at 1:500 dilution
Lane A: Rat brain tissue lysate
Lane B: Mouse brain tissue lysate

Lysates/proteins at 30 μ g per lane.
Secondary
Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.
Performed under reducing conditions.

Predicted band size: 98 kDa
Observed band size: 115 kDa