

ALK-1 / ACVRL1 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10066-R104

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human ALK-1 / ACVRL1 Protein (Catalog#10066-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human ALK-1 / ACVRL1 (rh ALK-1 / ACVRL1; Catalog#10066-H08H; NP_000011.2; Met1-Gln118).
Ig Type:	Rabbit IgG
Clone ID:	104
Specificity:	Human ALK-1 / ACVRL1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA(Cap)
RECOMMENDED CONCENTRATION	
Sandwich ELISA (Capture Ab)	ELISA(Cap): 1:250-1:2000 This antibody will detect Human ALK-1 / ACVRL1 in ELISA pair set (Catalog: # SEK10066). In a sandwich ELISA, it can be used as capture antibody when paired with (Catalog: # 10066-R024).

Please Note: Optimal concentrations/dilutions should be determined by the end user.