

AGRP Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10070-MM03

GENERAL INFORMATION

Immunogen:	Recombinant Human AGRP Protein (Catalog#10070-H01H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human AGRP (rh AGRP; Catalog#10070-H01H; NP_001129.1; Ser85-Thr132). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG2a
Clone ID:	03
Specificity:	Human AGRP
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000
-------	----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.