

Dppa2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 100950-T10

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to the center region of the Human Dppa2
Preparation	Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human Dppa2, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human Dppa2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	PESCRG1,ECAT15-2,CT100

APPLICATIONS

Applications:	IHC-P, ICC/IF
----------------------	---------------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:1000-1:5000
ICC/IF	ICC/IF: 1:300-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Dppa2 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist

Catalog Number: 100950-T10



Immunohistochemical staining of human Dppa2 in human brain with rabbit polyclonal antibody (1:2000, formalin-fixed paraffin embedded sections).



Immunofluorescence staining of Dppa2 in hESC-H9 cells. Cells were fixed with 4% PFA, permeabilized with 0.5% Triton X-100 in PBS, blocked with 10% serum, and incubated with rabbit anti-Human Dppa2 polyclonal antibody (dilution ratio 1:1000) at 4°C overnight. Then cells were stained with the Alexa Fluor 488-conjugated Goat Anti-rabbit IgG secondary antibody (green). Positive staining was localized to nucleus.