

CD34 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10103-MM17

GENERAL INFORMATION	
Immunogen:	Full length human CD34 expression plasmid
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with full length human CD34 expression plasmid. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	17
Specificity:	Human CD34
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	CD34
APPLICATIONS	
Applications:	IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:30-1:100


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD34 Antibody, Mouse MAb




Sino Biological
Biological Solution Specialist


Catalog Number: 10103-MM17




Immunohistochemical staining of human CD34 in human kidney with mouse monoclonal antibody at 1:60 dilution, formalin-fixed paraffin embedded sections.




Immunohistochemical staining of human CD34 in human corpus uteri with mouse monoclonal antibody at 1:60 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human CD34 in human stomach with mouse monoclonal antibody at 1:60 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human CD34 in human liver with mouse monoclonal antibody at 1:60 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human CD34 in human umbilical stalk with mouse monoclonal antibody at 1:60 dilution, formalin-fixed paraffin embedded sections.