

IL-1R1 / CD121a Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10126-R024

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human IL-1R1 / CD121a Protein (Catalog#10126-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human IL-1R1 / CD121a (rh IL-1R1 / CD121a; Catalog#10126-H08H; NP_000868.1; Met1-Thr332).
Ig Type:	Rabbit IgG
Clone ID:	024
Specificity:	Human IL-1R1 / CD121a
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human IL-1R1 / CD121a.

Please Note: Optimal concentrations/dilutions should be determined by the end user.