

Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) Matrix protein VP40 Antibody, Rabbit PAb, Antigen Affinity Purified
Catalog Number: 101298-T32



Sino Biological
Biological Solution Specialist


GENERAL INFORMATION	
Immunogen:	A synthetic peptide corresponding to the N-terminus of the Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) Matrix protein VP40
Preparation	Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) Matrix protein VP40 , and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) Matrix protein VP40
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	WB
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
Western Blot	WB: 1:1000-1:5000

Please Note: Optimal concentrations/dilutions should be determined by the end user.

**Ebola virus EBOV (subtype Zaire, strain
H.sapiens-wt/GIN/2014/Kissidougou-C15)
Matrix protein VP40 Antibody, Rabbit PAb,
Antigen Affinity Purified**
Catalog Number: 101298-T32



Sino Biological
Biological Solution Specialist



Anti-ZaireEBOV-GIN-AHX24648-VP40 rabbit
polyclonal antibody at 1:1000 dilution
Sample:ZaireEBOV-GIN-AHX24648-VP40
Recombinant Protein
Lane A: 20ng

Secondary
Dylight 800-labeled Antibody To Rabbit IgG
(H+L) at 1:5000 dilution.

Developed using the Odyssey technique.
Performed under reducing conditions.