

IL1B/IL-1B/IL-1 beta Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10139-MM39

GENERAL INFORMATION

Immunogen:	Recombinant Human IL1B/IL-1B/IL-1 beta Protein (Catalog#10139-HNAE)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human IL1B/IL-1B/IL-1 beta (rh IL1B/IL-1B/IL-1 beta; Catalog#10139-HNAE; NP_000567.1; Ala117-Ser269). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	39
Specificity:	Human IL1B/IL-1B/IL-1 beta
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	IHC-P
---------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:100-1:500
-------	--------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

IL1B/IL-1B/IL-1 beta Antibody, Mouse MAb




Sino Biological
Biological Solution Specialist

Catalog Number: 10139-MM39



Immunochemical staining of human IL1b in human lung with mouse monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.



Immunochemical staining of human IL1b in human liver with mouse monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.