

Endoglin / CD105 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10149-R103

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Human Endoglin / CD105 protein (Catalog#10149-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human Endoglin / CD105 (rh Endoglin / CD105 ; Catalog#10149-H08H; NP_001108225.1; Met1-Gly586).
Ig Type:	Rabbit IgG
Clone ID:	103
Specificity:	Human Endoglin / CD105
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	END,HHT1,ORW1

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION

Flow Cytometry	FCM: 1:25-1:100
-----------------------	-----------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.


Endoglin / CD105 Antibody, Rabbit MAb

Catalog Number: 10149-R103



Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.



Flow cytometric analysis of Human CD105 expression in HeLa cells.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.