

Endoglin/CD105 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10149-T24

GENERAL INFORMATION

Immunogen:	Recombinant Human Endoglin/CD105 Protein (Catalog#10149-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human Endoglin/CD105 (rh Endoglin/CD105; Catalog#10149-H08H; NP_001108225.1; Met1-Gly586). Endoglin/CD105 specific IgG was purified by Human Endoglin/CD105 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human Endoglin/CD105
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA, IHC-P
---------------	--------------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:1000-1:5000
ELISA	ELISA: 1:5000-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Endoglin/CD105 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist

Catalog Number: 10149-T24



Immunohistochemical staining of human CD105 in human kidney with rabbit polyclonal antibody at 1:2000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human CD105 in human placenta with rabbit polyclonal antibody at 1:2000 dilution, formalin-fixed paraffin embedded sections.