

Complement Component C2 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10154-MM06

GENERAL INFORMATION

Immunogen:	Recombinant Human Complement Component C2 Protein (Catalog#10154-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human Complement Component C2 (rh Complement Component C2; Catalog#10154-H08H; NP_000054.2; Met1-Leu752). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	06
Specificity:	Human Complement Component C2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA(Det)
----------------------	------------

RECOMMENDED CONCENTRATION

Sandwich ELISA (Detection Ab)	ELISA(Det): 1:1000-1:10000 This antibody will detect Human Complement Component C2 in ELISA pair set (Catalog: # SEK10154). In a sandwich ELISA, it can be used as detection antibody when paired with (Catalog: # 10154-R014).
--	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.