

CNTN4 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10178-T24

GENERAL INFORMATION

Immunogen:	Recombinant Human CNTN4 Protein (Catalog#10178-H08B)
Preparation	Produced in rabbits immunized with purified, recombinant Human CNTN4 (rh CNTN4; Catalog#10178-H08B; Q8IWV2-1; Met1-Ser1000). CNTN4 specific IgG was purified by Human CNTN4 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CNTN4
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA, IHC-P
----------------------	--------------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:50-1:200
ELISA	ELISA: 1:5000-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CNTN4 Antibody, Rabbit PAb, Antigen Affinity Purified


Catalog Number: 10178-T24



Sino Biological
Biological Solution Specialist



Immunohistochemical staining of human CNTN4 in human kidney with rabbit polyclonal antibody at 1:100 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human CNTN4 in human brain with rabbit polyclonal antibody at 1:100 dilution, formalin-fixed paraffin embedded sections.