

CD147 / EMMPRIN / Basigin Antibody (HRP), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10186-R118-H

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human CD147 / EMMPRIN / Basigin Protein (Catalog#10186-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CD147 / EMMPRIN / Basigin (rh CD147 / EMMPRIN / Basigin; Catalog#10186-H08H; NP_940991.1; Met1-His205), and conjugated with horseradish-peroxidase (HRP).
Ig Type:	Rabbit IgG
Clone ID:	118
Specificity:	Human CD147 / EMMPRIN / Basigin
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	BSG
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 0.1-1 µg/ml

Please Note: Optimal concentrations/dilutions should be determined by the end user.