

CD300f Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 101913-T32

GENERAL INFORMATION	
Immunogen:	A synthetic peptide corresponding to the C-terminus of the Rat CD300f
Preparation	Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Rat CD300f, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Rat Mouse (Species predicted to react based on 100% sequence homology)
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	WB
RECOMMENDED CONCENTRATION	
Western Blot	WB: 1:500-1:2000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD300f Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 101913-T32



Anti-Cd300f rabbit polyclonal antibody at
1:500 dilution

Lane A: Rat heart tissue lysate

Lane B: Rat liver tissue lysate

Lane C: Rat brain tissue lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000
dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 36 kDa

Observed band size: 36 kDa