

E-Cadherin / CDH1 / E-cad / CD324 Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10204-MM08-P

GENERAL INFORMATION

Immunogen:	Recombinant Human Cadherin-1 / CDH1 / CD324 / E-cadherin protein (Catalog#10204-H08H)
Reagents:	PE-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human Cadherin-1 / CDH1 / CD324 / E-cadherin (rh Cadherin-1 / CDH1 / CD324 / E-cadherin; Catalog#10204-H08H; NP_004351.1; Met 1-Ile 707) and conjugated with PE under optimum conditions, the unreacted PE was removed.
Ig Type:	Mouse IgG2a
Clone ID:	08
Specificity:	Human Cadherin-1 / CDH1 / CD324 / E-cadherin
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.09% sodium azide
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Alternative Names:	CDH1

APPLICATIONS

Applications:	FCM
----------------------	-----


Please Note: Optimal concentrations/dilutions should be determined by the end user.

E-Cadherin / CDH1 / E-cad / CD324 Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10204-MM08-P



Profile of anti-E-cad (CD324) reactivity on
MCF-7 cells analyzed by flow cytometry.

Flow cytometry was performed on a BD
FACSCalibur flow cytometry system.
Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for
technical protocols.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>