



Catalog Number: 10208-R107

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Human Layilin protein (Catalog#10208-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human Layilin (rh Layilin; Catalog#10208-H08H; NP_849156.1; Met 1-Glu 220).
Ig Type:	Rabbit IgG
Clone ID:	107
Specificity:	Human Layilin
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human Layilin.
--------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.