

**CD95 / APO-1 / TNFRSF6 / FAS
Antibody, Rabbit PAb, Antigen
Affinity Purified**

Catalog Number: 10217-T16



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION	
Immunogen:	Recombinant Human CD95 / APO-1 / TNFRSF6 / FAS Protein (Catalog#10217-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human CD95 / APO-1 / TNFRSF6 / FAS (rh CD95 / APO-1 / TNFRSF6 / FAS; Catalog#10217-H08H; NP_000034.1; Met1-Asn173). CD95 / APO-1 / TNFRSF6 / FAS specific IgG was purified by Human CD95 / APO-1 / TNFRSF6 / FAS affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CD95 / APO-1 / TNFRSF6 / FAS
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human CD95 / APO-1 / TNFRSF6 / FAS.

Please Note: Optimal concentrations/dilutions should be determined by the end user.