

G-CSFR / CD114 / CSF3R Antibody (APC), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10218-MM06-A

GENERAL INFORMATION

Immunogen:	Recombinant Human G-CSFR / CD114 / CSF3R Protein (Catalog#10218-HCCH)
Reagents:	APC-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human G-CSFR / CD114 / CSF3R (rh G-CSFR / CD114 / CSF3R; Catalog#10218-HCCH; NP_000751.1; Met1-Pro621) and conjugated with APC under optimum conditions, the unreacted APC was removed.
Ig Type:	Mouse IgG1
Clone ID:	06
Specificity:	Human G-CSFR / CD114 / CSF3R
Concentration:	10 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.09% sodium azide
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Alternative Names:	CSF3R

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

G-CSFR / CD114 / CSF3R Antibody (APC), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10218-MM06-A



Flow cytometric analysis of anti-human GCSFR (CD114) on human whole blood granulocytes. The fluorescence histograms were derived from events with the forward and side light-scatter characteristics of viable granulocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.