

# Complement C1s Antibody, Rabbit MAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 10220-R006

EliteRmab® is a registered trademark of Sino Biological Inc.

| GENERAL INFORMATION       |                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen:                | Recombinant Human Complement C1s protein (Catalog#10220-H08H)                                                                                                                                                                                      |
| Preparation               | This antibody was obtained from a rabbit immunized with purified, recombinant Human Complement C1s (rh Complement C1s; Catalog#10220-H08H; NP_001725.1; Met 1-Asp 688).                                                                            |
| Ig Type:                  | Rabbit IgG                                                                                                                                                                                                                                         |
| Clone ID:                 | 006                                                                                                                                                                                                                                                |
| Specificity:              | Human Complement C1s                                                                                                                                                                                                                               |
| Formulation:              | 0.2 µm filtered solution in PBS                                                                                                                                                                                                                    |
| Storage:                  | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles. |
| APPLICATIONS              |                                                                                                                                                                                                                                                    |
| Applications:             | ELISA                                                                                                                                                                                                                                              |
| RECOMMENDED CONCENTRATION |                                                                                                                                                                                                                                                    |
| ELISA                     | ELISA: 1:5000-1:10000<br>This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human C1S.                                                                                                                  |

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**