

Leptin / LEP Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10221-RP02

GENERAL INFORMATION	
Immunogen:	Recombinant Human Leptin protein (Catalog#10221-HNAE)
Preparation	Produced in rabbits immunized with purified, recombinant Human Leptin (rh Leptin; Catalog#10221-HNAE; NP_000221.1; Val 22-Cys 167). Leptin specific IgG was purified by human Leptin affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human Leptin / LEP
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA,IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:500-1:2000
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human Leptin.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Leptin / LEP Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 10221-RP02



Sino Biological
Biological Solution Specialist



Immunohistochemical staining of human LEP in human stomach with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).