

KRT10 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 102250-T08

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to the C-terminus of the Human KRT10
Preparation	Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human KRT10, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse, Rat
	Human, Mouse, Rat
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	IHC-P
----------------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:100-1:500
--------------	--------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

KRT10 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist


Catalog Number: 102250-T08




Immunohistochemical staining of KRT10 in human skin with rabbit polyclonal antibody at 1:300 dilution, formalin-fixed paraffin embedded sections.




Immunohistochemical staining of KRT10 in human breast carcinoma with rabbit polyclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of KRT10 in human lung cancer with rabbit polyclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of KRT10 in mouse skin with rabbit polyclonal antibody at 1:300 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of KRT10 in rat skin with rabbit polyclonal antibody at 1:300 dilution, formalin-fixed paraffin embedded sections.