

Cadherin-4 / CDH4 / R-CAD Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10230-RP01

GENERAL INFORMATION

Immunogen:	Recombinant Human CDH4 protein (Catalog#10230-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human CDH4 extracellular domain (rh CDH4; Catalog#10230-H08H; NP_001785.2; Met 1-Ala 734). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human Cadherin-4 / CDH4
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human CDH4.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.