

EphB4 Antibody (APC), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10235-R001-A

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Human EphB4 protein (Catalog#10235-H08H)
Reagents:	APC-conjugated Rabbit monoclonal antibody
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human EphB4 (rh EphB4; Catalog#10235-H08H; NP_004435.3; Met 1-Ala 539) and conjugated with APC under optimum conditions, the unreacted APC was removed.
Ig Type:	Rabbit IgG
Clone ID:	001
Specificity:	Human EphB4
Concentration:	10 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2℃-8℃ for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !
Alternative Names:	EphB4

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.

EphB4 Antibody (APC), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10235-R001-A

EliteRmab® is a registered trademark of Sino Biological Inc.

Flow cytometric analysis of anti-human EphB4 on MCF-7 cells. MCF-7 cells were detached using 1X trypsin, washed, then stained with APC Rabbit anti-human EphB4.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.