

METAP1 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10241-T48

GENERAL INFORMATION

Immunogen:	Recombinant Human METAP1 Protein (Catalog#10241-H08B)
Preparation	Produced in rabbits immunized with purified, recombinant Human METAP1 (rh METAP1; Catalog#10241-H08B; NP_055958.2; His52-Phe386). METAP1 specific IgG was purified by Human METAP1 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human METAP1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB, ELISA
----------------------	-----------

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

METAP1 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10241-T48



Anti-METAP1 rabbit polyclonal antibody at
1:500 dilution
Lane A: A549 Whole Cell Lysate

Lysates/proteins at 30 µg per lane.
Secondary
Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000
dilution.

Developed using the ECL technique.
Performed under reducing conditions.

Predicted band size: 43 kDa
Observed band size: 43 kDa