

GH2/Growth hormone 2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 102422-T40

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to the N-terminus of the Human GH2/Growth hormone 2
Preparation	Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human GH2/Growth hormone 2, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human GH2/Growth hormone 2
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB, IHC-P
----------------------	-----------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:100-1:500
Western Blot	WB: 1:500-1:2000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

GH2/Growth hormone 2 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist


Catalog Number: 102422-T40



Immunohistochemical staining of human GH2 in human placenta with rabbit polyclonal antibody at 1:300 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human GH2 in human colon carcinoma with rabbit polyclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.



Anti-GH2 rabbit polyclonal antibody at 1:500 dilution

Lane A: H1299 Whole Cell Lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 25 kDa

Observed band size: 25 kDa