

Fas Ligand/FASLG Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10244-T16

GENERAL INFORMATION

Immunogen:	Recombinant Human Fas Ligand/FASLG Protein (Catalog#10244-H07Y)
Preparation	Produced in rabbits immunized with purified, recombinant Human Fas Ligand/FASLG (rh Fas Ligand/FASLG; Catalog#10244-H07Y; NP_000630.1; Pro134-Leu281). Fas Ligand/FASLG specific IgG was purified by Human Fas Ligand/FASLG affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human Fas Ligand/FASLG
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human Fas Ligand/FASLG.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.